



## Kleenpak Presto sterile connectors

Genderless sterile connections

### Overview

**The Kleenpak Presto sterile connector is a genderless connection that allows for the permanent sterile connection of two fluid streams for a large range of biopharmaceutical applications.**

With its intuitive operation, sterile connections can be carried out in a simple three-step operation, even in an unclassified environment, without compromising the sterility of the fluid paths.

Key design features such as the anti-actuation tabs and tamper-resistant protective cap give end users more confidence in operation and security of supply.

The Kleenpak Presto sterile connector brings enhanced levels of quality assurance with 100% inspection at point of manufacture that allows device traceability, as each device has its own serial number. The manufacturing process of Kleenpak Presto sterile connectors brings an enhanced level of quality assurance, with the use of an automatic vision system that ensures the absence of defects in the membrane and membrane welding.

Available in a variety of sizes – 6.35 mm ( $\frac{1}{4}$  in.), 9.53 mm ( $\frac{3}{8}$  in.), 12.7 mm ( $\frac{1}{2}$  in.), 15.8 mm ( $\frac{5}{8}$  in.), 19 mm ( $\frac{3}{4}$  in.) hose barb, and 12.7 mm ( $\frac{1}{2}$  in.) sanitary connection – the Kleenpak Presto sterile connector can be used in upstream processing, downstream processing, and formulation and filling. Made from Bisphenol-A (BPA) free Polyethersulfone (PES) the Kleenpak Presto sterile connector is compatible with a wide range of process fluids and solvents.

### Applications

Kleenpak Presto Sterile Connector

Media preparation and transfer

Buffer preparation and transfer

Transfer of inoculum to bioreactor

Sampling during fermentation / cell culture

Bioreactor harvest

Sterile fluid transfer between unit operations

Bulk handling of sterile material in non-classified environments

Probe insertion into bioreactors, mixers and 3D biocontainers

Sterile filtration manifolds

Hybrid stainless steel and single-use system connection

Connection of bulk sterile material to filling machine

Sterile waste removal from process streams

## Quality Standards

Manufactured under a quality management system certified to ISO9001

Manufactured in a clean room Class 7 in operation

Supplied with a certificate of test confirming the quality standards and quality control tests performed by Pall

Each connector is individually marked with batch number and serial number

Batch release criteria:

USP 85 - Endotoxin

USP 788 - Particulate test

Device release criteria:

100% inspection through vision system for absence of membrane and weld defects.

The fluid path materials of construction have been tested and meet the regulatory requirements of:

USP 88 - Biological reactivity in vivo for Class VI 121°C plastics

USP 87 - Biological reactivity in vitro

USP 661 - Physicochemical tests

The fluid path materials of construction do not contain substances derived from animal products (i.e. BSE/TSE risk free)

## Validation Tests

Mechanical tests

Leak (closure integrity); burst; creep rupture and pressure hold test; tensile strength

Functional tests

Extreme temperature pressure leak test strength and performance testing

Water flow characteristics

Bacterial challenge ("soiling") test using *Brevundimonas diminuta*

Autoclave and gamma resistance

Shelf life studies

Extractables testing

## Product specifications

### ½ in. hose barb connection, box of 50

Catalog code	PSC1G06
Price	2,961.00 USD
Connection	½ in hose barb
Connector Body Material	Polyethersulfone (PES) (in fluid path)
Connector Format	Genderless
Internal Diameter (Metric)	14 mm
Max. Gamma Irradiation Dose	50 kGy
Max. operating pressure	3 barg up to 90 days; 4 barg up to 2 days

<b>Max. Operating Temperature</b>	60°
<b>Min. Operating Temperature</b>	2°
<b>O-ring material</b>	Platinum-cured silicone (in fluid path)
<b>Sterilization Conditions (Autoclave)</b>	75 minutes at 130°C