

## **DESCRIPTION:**

Saint-Gobain VueLife® "AC" Series Bags utilize a treated form of fluorinated ethylene propylene (FEP) film. This proprietary surface treatment results in the VueLife "AC" Bag interior having a higher surface energy than the VueLife® "C" Bag. The modified surface of the VueLife "AC" Bag is meant to better promote cell and protein adhesion to the surface of the bag, a feature that can be leveraged for transduction and transfection applications. VueLife bags are made from the highest quality USP Class VI FEP material and provides permeability to oxygen and carbon dioxide whilst remaining impermeable to water for improved culture and expansion.

Each "AC" series bag typically features a needle-less injection site and a "Y" connector with PVC tubing leading to a female luer and a heat-sealed sterile docking tube to greatly reduce the risk of culture contamination while allowing easy access for the introduction of new material. Other types of ports, tubing and connections are available and can be customized to the needs of the user.

### **APPLICATION:**

Adhesion properties allow for the coating of transduction enhancers.

#### FEATURES AND BENEFITS:

- <u>Full Extractables Report per BioPhorum</u> Operations Group (BPOG) Protocol
- <u>Technical Dossier Available</u>
- High purity
- Closed system

- Extreme Operating Temperature Range
- Chemical and Biologically Inert
- Gas Permeable
- Transparent


VueLife <sup>®</sup> "AC" Bag Specifications and Technical Information						
Component	Material	Length/Size				
Bag Fluid Contact Layer	5 Mil FEP	See VueLife <sup>®</sup> "C" Bag Size Specifications Chart				
Female Luer Barb Connector	PVDF	Female Luer w/ Barb for 0.125" (3.2mm) ID Tubing				
Small Pinch Clamp	PP	For up to .25" (6.3mm) OD Tubing				
3-Way Solvent Bond Y	PVC	0.161" (4.1mm) OD socket joints				
Male Luer Barb Connector	PVDF	Male Luer w/ Barb for 0.125" (3.2mm) ID Tubing				
FEP Retaining Ring	FEP	Ring made from .25" (6.3mm) OD x 0.1875" (4.8mm) ID Tubing				
Needleless Injection Site	Polycarbonate/Silicone	0.142" (3.6mm) ID				
Female Luer Port	FEP	Female Luer w/ 0.125" (3.2mm) ID				
Tubing	PVC	0.158" (4mm) OD x 0.118" (3mm) ID				
Spike Port with Port Pocket*	FEP	0.204" (5.2mm) ID				

\*750-C1 only

Measurements are for reference only

VueLife <sup>®</sup> "AC" Bag Size Specifications						
Product	Recommended Process Volume	*Approximate Maximum Fill Volume	Inner Bag Length in (cm)	Inner Bag Width in (cm)	Approximate Surface Area (x2 Bag Face)	
32-AC	32 mL	45 mL	2.1" (5.3)	3.2" (8.1)	87 cm <sup>2</sup>	
72-AC	72 mL	185 mL	4.4" (11.2)	3.4" (8.6)	193 cm <sup>2</sup>	
118-AC	118 mL	285 mL	4.3" (10.9)	5″ (12.7)	277 cm <sup>2</sup>	
119-AC	119 mL	290 mL	4.1′ (10.4)	5.4" (13.7)	286 cm <sup>2</sup>	
197-AC	197 mL	505 mL	7.2" (18.3)	5.4" (13.7)	502 cm <sup>2</sup>	
290-AC	290 mL	1000 mL	9.9" (25.1)	5.2" (13.1)	658 cm <sup>2</sup>	
750-AC	750 mL	2500 mL	15.4" (39.2)	8.3" (21.0)	1643 cm <sup>2</sup>	

\*The maximum fill volume is a calculated value for the theoretical maximum volume based on the bag dimensions. There is no data supporting usage at this volume. The recommend process volume for cell culture is provided to achieve the optimal 1 cm height when the bag is lying flat. Using a volume above the recommended process volume has not been studied.



### **ORDERING INFORMATION:**

Normally not stocked, built to order Sold as 10 individually packaged sterile bags per pack.

# Saint-Gobain Life Sciences:

50 W Watkins Mill Rd Gaithersburg, MD 20878

Customers should validate this product for their specific application

IMPORTANT: It is the user's responsibility to ensure the suitability and safety of Saint-Gobain Life Sciences products for all intended uses and that the materials to be used comply with all applicable medical regulatory requirements. Saint-Gobain Life Sciences assumes no responsibility for any product failures that occur due to misuse of the materials it provides arising out of the design, fabrication or application of the products into which the materials are incorporated.

WARRANTY: For a period of 12 months from the date of first sale, Saint-Gobain Life Sciences warrants this product to be free of defects in materials and workmanship. Our only obligation will be to replace any portion proving defective, or at our option, to refund the purchase price there of.

SAINT-GOBAIN LIFE SCIENCES DISCLAIMS ALL IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.