Item #



Job _____

STARSYS DIFFICULT AIRWAY CART

Areas such as the Emergency Department have varying case loads and requirements on a daily basis. In addition, they must mobilize or reset to address dramatic changes in process and flow, as experienced in disaster response. Starsys provides ultimate flexibility to reconfigure the space at a moments notice.

- Holds Manual Jet Ventilator with Regulator
- Tube Wrap
- $-O_2$ Tank Holder
- Enclosed Bronchoscope Holder

Storage for:

- Preparation
- Laryngoscope equipment accessories
- Endotracheal tubes and Styletes
- Supraglottic devices
- Surgical Airway devices
- Complete Storage System Addresses applications facility-wide.
- Interchangeable Components Provides flexibility to address changing needs.
- **Modular Design** Provides flexibility to change layout and/or location.
- **Polymer Construction** Maintains function and aesthetic over time.

Metro incorporates several elements in its product design to support a facility's cleaning protocols:

- Microban[®] antimicrobial protection helps prevent the growth of stain and odor causing bacteria on the product. Look for the "red check" symbol for this added protection.
- Advanced polymer and other proprietary finishes provide corrosion resistance
- Smooth rounded corners to allow for easier cleaning





Passive locking with plastic seals provides a level of visual security that indicates potential supply discrepancy without compromising access.

• Saves Time — Provides instant inventory check. Secures each drawer individually. Inventory only those drawers that have a broken seal.



- Increases Productivity Removable totes reduce down time by allowing instant replenishment with pre-stocked totes.
- **Improves Ergonomics** Drawers are easily reconfigured without tools allowing the most frequently used drawers to be relocated to the top of the cart.

*MICROBAN® and the MICROBAN® symbol are registered trademarks of the Microban Products Company, Huntersville, NC.




InterMetro Industries Corporation

North Washington Street Wilkes-Barre, PA 18705 www.metro.com Job











Right side of cart contents:

Harness to hold O₂ tank Hooks to hold jet ventilator

Left side of cart contents: Bougie stylete (coude) Airway exchanger (8fr and 14fr) Fiberoptic Lighted stylete (pediatric and adult size)

- regulator Hooks to hold jet vent catheters Enclosed Bronchoscope
- Holder

Dimensions Model No.

SXRSDIFAIR

 Overall Dimensions (Including Accessories) (in.) (mm)

 31 L. x 24.3 W. x 57.8 H.
 787 L. x 617W. x 1468H.





Drawer 1: Contents -**Preparation Drawer** (Working top left corner, down) Berman intubating oral airways (for intubating bronchoscope) Oral airways End-tidal CO₂ disposable detector Betadine solution Lidocaine ointment (2% Xylocaine) Hurricaine/Benzocaine topical anesthetic spray Oxymethalozone nasal spray Nasal airways Suction devices (i.e. suction tubing, yankauer, frazier suction and suction catheter) Bite blocks



Equipment/Accessories Assorted syringes Batteries for laryngoscope handles Assorted laryngoscope blades Disposable AirTraq device Laryngoscope handles



Drawer 3: Endotraceal Tubes and Styletes Disposable LMAs (Sizes 3, 4 and 5) LMA ET tubes (Sizes 6, 6.5, 7, 7.5 and 8) Stabilizer rods Styletes (adult and pediatric sizes)



Drawer 4: Supraglottic Device King airways (Size 3, 4 and 5)



Drawer 5: Surgical Airway Devices Melker Cricothyrotomy kits Surgical airway kits Tracheostomies, with cannulas (assorted sizes) Transtracheal Jet Ventilator regulator Transtracheal Jet Ventilation catheter (pediatric and adult sizes)

All Metro Catalog Sheets are available on our Web Site: www.metro.com



InterMetro Industries Corporation

North Washington Street, Wilkes-Barre, PA 18705 Product Information. U.S. and Canada: 1.800.992.1776 Outside U.S. and Canada: www.metro.com/contactus

L03-280 Printed in U.S.A. 9/10 Information and specifications are subject to change without notice. Please confirm at time of order.

Copyright © 2020 InterMetro Industries Corp