

TECHNIGLOVE DATA SHEET

TechNiGlove
♦ INTERNATIONAL ♦

RV400 & CP400

Low Extractable Powder Free Nitrile Examination Glove
Does not contain Thiuram, Carbamate and/or Thiazole

PHYSICAL PROPERTIES

Material	100% Nitrile
Style	Non-Sterile Ambidextrous Powder Free
Length	9.5" (240mm)
Thickness	4 mil (.004")
Grip Surface	Textured Fingertips
Cuff	Beaded
Color	The New Blue

	Before Aging	After Aging
Tensile Strength (MPa)	≥17	≥16
Elongation (%) Specification	650 min.	600 min.
Elongation (%) Typical Value	700	650
Force at Break (N)	≥6	≥6

ESD Specification @ RH45-55%, Temp: 20-26C
(Ref: ANSI/ESD- STM11.11/ANSI/ESD SPI5.1)

Surface Resistivity (Ohm/sq)	≤ / = 10 ¹¹
Decay Time (Second)	≤ / = 2.0
Tribo Charge (Volt)	≤ / = 50



**RiVAL™ – The New Blue
Accelerator Free Nitrile**

PACKAGING INFORMATION

9.5" Standard Boxed, Nitrile Gloves

RV400	X-SMALL
RV401	SMALL
RV402	MEDIUM
RV403	LARGE
RV404	X-LARGE
RV405	XX-LARGE

200 Gloves/Box - 2,000 Gloves/Case

9.5" Clean-Packed, Nitrile Gloves

CP400	X-SMALL
CP401	SMALL
CP402	MEDIUM
CP403	LARGE
CP404	X-LARGE
CP405	XX-LARGE

200 Gloves/Bag - 2,000 Gloves/Case

TechNiGlove
♦ INTERNATIONAL ♦

3750 Pierce Steet, Riverside, CA 92503 USA
Toll Free (877) 648-7453 Phone (951) 582-0890 Fax (951) 582-0990
www.TechNiGlove.com