

TECHNICAL DATA SHEET

TechNiGlove
♦ INTERNATIONAL ♦

RVCR700 & RVCR1500

Rival CR's new soft modulus nitrile delivers the fit and feel of latex for the ultimate stretch and comfort. Suitable for cleanroom, pharmaceutical, biotech and medical device manufacturing.

PHYSICAL PROPERTIES

Style	Non-Sterile Ambidextrous
Length	RVCR700 9.5" (240mm) RVCR1500 12" (300mm)
Thickness	4 mil (.004")
Grip Surface	Textured Fingertips
Cuff	Beaded
Color	White
Tensile Strength	> 14 MPa
Elongation at Break	> 650



**RiVAL™ – The New Blue
Class 100 Accelerator Free Nitrile**

TECHNICAL PROPERTIES

Particle Levels	<1800 total particles/cm ² IEST- RP-CC005.4 Method
Surface Resistance in Use	<1x10 ⁹ ohm/square per ANSI/ESD SP15.1

Extractables

TNVR	< 5.0 ug/cm ²
Fluoride (F ⁻)	< 0.50 ug/cm ²
Chloride (Cl ⁻)	< 1.5 ug/cm ²
Bromide (Br ⁻)	< 0.05 ug/cm ²
Nitrate (NO ₃ ⁻)	< 1.00 ug/cm ²
Phosphate (PO ₄ ³⁻)	< 0.05 ug/cm ²
Nitrite	< 0.50 ug/cm ²
Sulphate (SO ₄ ²⁻)	< 1.00 ug/cm ²

PACKAGING INFORMATION

9.5" Class 100 White Soft Modulus

		Case Wt
RVCR700	X-SMALL	10 lbs
RVCR701	SMALL	11 lbs
RVCR702	MEDIUM	11 lbs
RVCR703	LARGE	12 lbs
RVCR704	X-LARGE	12 lbs
RVCR705	XX-LARGE	13 lbs

100 Gloves/Poly Bag, 10 Poly Bags/Case = 1,000 Gloves/Case
Case Size: 10.5"x13"x12"

12" Class 100 White Soft Modulus

		Case Wt
RVCR1500	X-SMALL	12 lbs
RVCR1501	SMALL	12 lbs
RVCR1502	MEDIUM	12 lbs
RVCR1503	LARGE	13 lbs
RVCR1504	X-LARGE	14 lbs
RVCR1505	XX-LARGE	15 lbs

100 Gloves/Poly Bag, 10 Poly Bags/Case = 1,000 Gloves/Case
Case Size: 10.5"x13"x12"

TechNiGlove
♦ INTERNATIONAL ♦

3750 Pierce Steet, Riverside, CA 92503 USA
Toll Free (877) 648-7453 Phone (951) 582-0890 Fax (951) 582-0990
www.TechNiGlove.com